



Sani-bag®

Description:

Sterile kit to secure the exit-site and placement catheters. The kit is composed by a pocket, plasters, dressings, and, on request, by Combi caps.

Functions:

particular pocket welded with outside material, soft and absorbent and inside waterproof, does not permit the escape of liquids and prevents by penetration from outside, making easy the patient's abluion. The placement of catheters is done easily through an entrance hole drilled in the Pocket. The practical Pocket closure allows total protection of enlightenment, isolating it completely from the external environment.

Package:

medical thermoformed tray with medical paper and polyethylene film. On each pan is visible after sterility indicator, part number, manufacturer name, lot number and expiry date, method of sterilization, symbology

Intended use:

for dressing and exit-site protection use.
 Non-invasive surgery device with a mechanical barrier function.

Warnings:

sterility is guaranteed until unopened. Store in a dry place, protected from light, at room temperature.

Disposal:

the product used in the health sector must be disposed in according to the regulations

Code	Product	Pcs/Pack	Pcs/trasport carton
500017/B	Sani-bag	1	300